genu3 & primaDOC

Apparecchiature isocinetiche monoarticolari e pluriarticolari
Single-joint and multi-joint isokinetic dynamometers

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The most effective, most compact, simplest and cheapest solution for the problems related to the evaluation and recovery of muscular strength and resistance through an isokinetic approach.

The genu3 and primaDOC apparatuses are designed for daily use. The accuracy of the measurements and the electronic feedback enhance the quality of the rehabilitation work and patient's motivation. The user-friendliness of these apparatuses cuts operator's time and thus increases the efficiency of the rehabilitation center.

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easytech for sport professionals

Roberto Baggio during a workout on the isokinetic **genu3** apparatus.

The isokinetic apparatuses are designed to permit the muscles to be exercised at a speed that remains constant through the entire range of motion of the joint. The resistance provided by the apparatuses is constantly proportional to the force applied by the patient. In this way, the involved muscles always develop the highest tension and maximum stimulation of the fibers.

These characteristics of the isokinetic apparatuses make them very versatile tools that can be advantageously utilized to evaluate the patients from the functional viewpoint, for rehabilitation and training.

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General characteristics

*The **genu3** and **primaDOC** apparatuses use an electronically controlled hydraulic brake to ensure best smoothness of the movement and full safety for the patient. In order to obtain a very compact design and reduce overall dimensions, the brake is included in the structure that supports the seat. The seat cover can be completely removed if necessary. The control console can be reached also from the isokinetic apparatus.*

***Genu3** and **primaDOC** are complete systems for rehabilitation, training and patient's evaluation. All the measured parameters are visualized in real time and can be stored for subsequent processing.*

*The software is developed in Windows™ environment, and is easy to use. It also helps the operator motivate the patient during the session thanks to the instantaneous feedback and an amusing video game. A special option enables the operator to connect with **easytech** for remote technical assistance.*

*Some software-controlled constraints ensure maximum internal and application safety (the joint range of motion can be limited through the adoption of "software stops" that can be placed at any possible angle (**genu3** is also fitted with mechanical stops).*

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Types of movement

- Isokinetic
In this case, the movement is at a constant speed, and the resistance provided by the apparatus changes as a function of the strength that the patient applies
- Inertial isotonic
The characteristic of the inertial, isotonic movement is the constant resistance provided by the apparatus. In this case, the strength that the patient applies changes as a function of the lever arm and inertia.
- Hydrodynamic
Hydrodynamic movement is obtained when the apparatus produces a resistance proportional to the square of the speed at which the movement is performed (viscous resistance) simulating the movement in the water.

Main parameters that can be measured

- Speed
- Torque
- Range of motion
- Work
- Power
- Index of resistance (fatigue)

The computer of the apparatus measures the parameters. The results are shown in numerical and graphical form, and are entered automatically in the pre-set medical report.

Models

genu3 is the single-joint system specifically designed to treat the knee; the special **genuV** bench (optional) permits exercises to be performed in Closed Kinetic Chain

primaDOC is the multi-joint system that lets you treat the following articulations:

- Knee (flexion/extension)
- Ankle (dorsal-plantar flexion)
- Shoulder (intra-extra rotation in abduction)
- Closed Kinetic Chain of the lower limb

Apparecchiature isocinetiche monoarticolari e pluriarticolari Single-joint and multi-joint isokinetic dynamometers

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Technical characteristics

	genu3	primaDOC
Configuration	Isokinetic unit: single-seat with integral oil brake. Seat cover can be fully removed Control console / processing unit: Pentium processor, or compatible, Windows, 15 in. screen, keyboard, mouse and color printer; connection to isokinetic unit type RS232 serial	Isokinetic unit: single-seat with integral oil brake. Seat cover can be fully removed Control console / processing unit: Pentium processor, or compatible, Windows, 15 in. Screen, keyboard, mouse and color printer; connection to isokinetic unit type RS232 serial Accessory basket: It can be stowed inside the cart/console
Movement control	Speed control range: 40°/sec – 400°/sec. Independent and separate adjustment for agonist – antagonist movement. Torque: measured over the entire joint angle, resolution 1 Nm Range of motion: 240°, internal resolution 0.75° (double values for shoulder angular device), measurement over the entire range Work: measurement in real time; maximum work, work per repetition, total work. Resistance index Power: real time measurement; maximum power; power per repetition, average power	
Operating configurations and options	Basic configuration: lower limb Option genuV: for exercises in Closed Kinetic Chain	Basic configuration: lower limb Extended configuration: lower limb and shoulder Shoulder angular device: it permits the change from the basic to the extended configuration
Overall dimensions and weights	Isokinetic unit: cm 110X110X100 h Control console: cm 70X70X115 h Minimum operating space: cm 200X270 Isokinetic unit: Kg 155 Control console: Kg 30 Option genuV: Weight: kg 36 • Overall dimensions: cm 132X35X135 h • Minimum operating space: genu3 + genuV, cm 200X520	Isokinetic unit: cm 70x125x125 h Control console: cm 60x57x92 h Minimum operating space: cm 150x200 Isokinetic unit: Kg 128 Control console: Kg 30
Electrical characteristics	Absorbed power: max 500VA (standby/operation). Protection fuses accessible from outside the apparatus Class and type according to standard CEI EN 60601: Class I - type B	
Installation requirements	Electrical power supply: Single-phase, 220/230 V – 50/60 Hz; available at 2 wall-mounted mains outlets Other supplies: Telephone line with cable and standard modular plug for analog telephone line (necessary for remote assistance option only) Environmental conditions: Ambient temperature from + 10° C to + 40° C. Relative humidity: 0% - 75% (without moisture)	
Standards	Apparatus conforming to directive CEE 93/42 (CE mark)	