


The logo for 'easytech' features the brand name in a bold, blue, sans-serif font. A stylized blue arc curves over the 'y' and under the 't', suggesting a medical or technological theme.

**easytech**

A photograph showing a therapist in a light grey polo shirt and white shorts examining a patient's leg. The therapist's hands are positioned to palpate the knee and lower thigh area. The patient's leg is extended horizontally across the frame.

Hydraulic and electric two sections treatment tables for medical examination.  
Practical and easy to set up, developed especially for examinations, treatments and  
massages.  
Designed to perfection to guarantee your patients maximum comfort during therapy.

## 2-section treatment tables



## 2-section treatment tables

Two sections treatment tables developed especially to take care of your patients.

Practical and easy to use, they are suitable for both examination and treatment.

**Hydraulic treatment tables** activate without any plugs, allowing the operator to adjust height by simply pressing the mechanical lever, placed underneath.

**Electric treatment tables** save your time and guarantee easy transportation. Although their use depends on electricity, they are able to offer your patients maximum comfort and pleasure.

**Those treatment tables are a class I medical devices** according to standard N.°13, section III, attachment VIII of EU regulation 2017/745.

### OVERALL SIZES:

Upholstery size: 1900 x 610 mm

Height regulation: 545 a 910 mm

Adjustment angle of the headboard: between 0° and 70°

### TECHNICAL FEATURES OF THE STRUCTURE:

- Welded steel pipe structure treated with anti-bacterial epoxy paint;
- Height regulation for the hydraulic model: hydraulic regulation by pressing on the pedal;
- Height regulation for the electric model: motore elettrico con pulsante da piede;
- Regulation of the back through a lever-operated gas spring on the back of the backrest;
- 4 50mm diameter polyurethane wheels with securing system;
- pedal for wheels retraction with retractable system;



Two sections hydraulic treatment table

- Leveling system for ground mounts;
- Optional paper holder suitable for 500mm paper rolls.

### UPHOLSTERY FEATURES :

- Upholstery divided in MDF with a 50 mm foam;
- Face opening for better patient comfort;
- Washable lining made of high strength nappa;
- Additional upholstery colors are available on request.

### SAFETY INSTRUCTION:

- The treatment table must be used only for the purpose it was intended for;
- The treatment table cannot be used for patients who weigh more than 180 kg;
- Attention must be paid to avoid flattening whenever the hydraulic system and its mobile parts are activated;
- The patient should seat exclusively in the specific seating area, avoiding this way any possible imbalance in the table and possible fall;
- The climbing of the patient onto the table must be carried out with the frame in its lowest position.

### STORAGE AND CONSERVATION:

- Every Hcaresol treatment table is packed in cardboard, should not be exposed to rain or placed in places with high humidity levels, to avoid any possible damage to the packaging or the product;
- The products should be stored in safe places that guarantee protection against falls or any other situation that could damage them;
- Before starting any type of cleaning operations consult our product instructions.



Two sections electric treatment table